# Journal of FisheriesSciences.com

E-ISSN 1307-234X

© 2017 www.fisheriessciences.com

**Research Article** 

### **Comparison of DNA-extraction Methods Suitable for PCR-based Applications to Identify Shrimp Species in Commercial Products**

Besbes Nadia and Saloua Sadok<sup>1</sup>

<sup>1</sup>Institut National des Sciences et Technologies de la Mer (INSTM), Laboratory of Blue Biotechnology and Aquatiques Bioproducts B3Aqua, Tunisia

Received: 02.10.2017 / Accepted: 21.11.2017 / Published online: 24.11.2017

Abstract: Seven different methods of DNA extraction were performed and compared in order to have the better DNA quality and quantity from specific seafood matrix such as various types of shrimp products. Thus, the nucleic acids extraction efficiency was evaluated according to DNA purity for PCR amplification. To evaluate the interaction of the different parameters, a multivariate regression, PLS approach was used to treat the results of PCR amplification of fresh and processed shrimp. The principal component (PCs) characterized the PCR amplification of shrimp (62% of all variance) with 59.0 and 3.0% for PC1 and PC2 respectively. Based on PC1 analysis, a positive correlation was obtained between the extraction procedures using PCI/SDS, DNeasy Tissue Kit (QIAGEN), TRITON and CTAB, and the results of PCR amplification. In addition to the successful optimization of DNA yield, the multivariate analysis revealed the different factors that affected the PCR amplification including the type of sample matrix. The comparison of the seven extractions methods, in this paper, has highlighted the efficacy of the SDS/PCI precipitation method for fresh and processed shrimp.

**Keywords:** Polymerase chain reaction; Processed shrimp; Mitochondrial DNA; Cytochrome b gene; Partial least square PLS1 regression

Sadok Saloua, Institut National des Sciences et Technologies de la Mer Port de Pêche La Goulette, La Goulette 2060, Tunisia, Tel: 00216 71 73 58 48; E-mail: salwa.sadok@instm.rnrt.tn

#### Introduction

The prevalence of mitochondrial DNA use for species differentiation and foodstuffs traceability is fairly well established (Chandrika et al., 2010). Although nuclear DNA is more powerful, mtDNA was more extensively exploited as it is more stable over time and conditions due to its circular configuration, presence in thousands of copies within mitochondria and is extractable from highly degraded source such as processed food (Brown et al., 1979, Wiesne et al., 1991, Prado et al., 2005, Kaniak-Golik et al., 2015).

Various different methods can be applied for DNA extraction from diverse sources of samples, but the recuperation of highquality of DNA represents a challenging step for subsequent applications such as electrophoresis, polymerase chain reaction (PCR) and sequencing. Difficulty in isolating intact and sufficient DNA from samples differs with the physico-biochemical natures of the tissue (Rumpho et al., 1994). For instance, polysaccharides, naturally present in tissue may not completely be removed using classical DNA extraction, causing inhibition of DNA amplification (Cawthorn et al., 2011; Besbes et al., 2015). Thus, DNA-based methods are highly dependent on the extraction and purification techniques (Lockley et al., 2002; Di Pinto et al., 2007) and selection of the appropriate method according to sample matrix type is crucial. For instance, for canned pelagic fish such as tune and sardines, sample pre-treatment including tissues defatting is necessary before DNA extraction (Bauer et al., 2003; Jérôme et al., 2003; Chapela et al., 2007).

During the last decades, various protocols have been described including the phenol-chloroform-based approaches (Besbes et al., 2011), the CTAB extraction method with modifications and the use of commercial kit (Jérôme et al., 2003, Jérôme et al., 2003, Olexová et al., 2004). For the specific shrimp matrix, seven DNA extractions were performed and were followed by purification procedures. For the comparative investigation, discussion was made on the basis of the best extraction, DNA purity and DNA suitability for amplification from fresh and processed shrimp. The aim of this work is to define the adequate, but relatively easy, reliable and cheap method to extract high quality and quantity of DNA from various types of shrimp products, for subsequent application including PCR and sequencing.

#### **Materials and Methods**

#### **Commercial food products**

Freshly caught shrimp (*Parapenaeus longirostris*, *Metapenaeus kerathurus*) along the Tunisian coast and refrigerated shrimp (*Metapenaeus monoceros*, *Plezionika heterocarpus*) bought from seafood plant; were taxonomically identified with respect to external characters. Three shrimp-type products were purchased from Tunisian supermarkets in frozen state (cooked, peeled precooked and breaded). Upon arrival in the laboratory, all samples were stored at -80°C until analysis. The extraction of mt-DNA from the flesh (cephalo-thorax part) was performed for all samples according to the different protocols (EX1- EX7). Careful attention of flesh sampling was given for breaded shrimp to avoid matrix contamination.

#### **DNA Extraction**

#### Kit fast DNA (EX1)

The DNA was extracted following the Kit -manufacturer's instructions. To the weighed subsample (0.1 g); 800  $\mu$ l Genomic Lyse buffer was added (to lyse the cells and solubilise the proteins, DNA and other cellular constituents) in the presence 0.1 mg (3.2 U) of Proteinase K (to digest the protein) and incubated at 65°C for 10 to 30 minutes. The DNA-containing solution was mixed with the genomic binding buffer (100  $\mu$ l), then eluated on a column which binds the DNA.

### Method using the phenol-chloroform-isoamyl alcohol: SDS/PCI (EX2)

The extraction of DNA was performed according to Besbes et al., 2011 with the following modifications. As shrimp muscle contains fewer lipids, a one step DNA extraction was performed. Sample (0.1 g) was homogenised with 500  $\mu$ L of the buffer A (50 mM Tris–HCl, 100 mM EDTA and 1% SDS with pH 8.0) and 0.1 mg (3.2 U) of proteinase K (Sigma Aldrich, France). The mixture was incubated at 55°C with continuous shaking during 4 h, cooled on ice and then centrifuged at 12,000g for 10 min. The supernatant was extracted with phenol–chloroform–isoamylic alcohol in a 25:24:1 ratio, then precipitated with 600  $\mu$ L of isopropanol at -20°C and mixed by inversion. Using ethanol (70%), the pellet was washed and dried, then resuspended in 50  $\mu$ L TE buffer.

#### Method using the SDS/PCI/NaCl (EX3)

In this case, the same procedure was applied to 0.1 g of sample using the SDS/PCI method (Hsieh et al., 2005), but with the addition of NaCl to the buffer A (final concentration of 0.2 M NaCl).

#### Method using Chelex (EX4)

The sample (0.1 g) was homogenised in 300  $\mu$ L of a Chelex water solution (5% -Chelex 100 resin Qiagen) mixed with 20  $\mu$ L of proteinase K (0.1 mg) and 30  $\mu$ L TE buffer; incubated at 56°C during 4 h and followed by a heating (100°C) for 15 min (Jérôme et al., 2003). The solution was centrifuged (12,000 g for 10 min), the DNA pellet for PCR amplification was stored at 4°C until use.

#### Method using CTAB/CI (EX5)

The shrimp sample (0.1 g) was grinded in 500  $\mu$ L pH 8-CTAB buffer (100 mM Tris–HCl, 20 mM EDTA, 2% hexadecyl– trimethyl ammonium bromide, and 1.4 M NaCl), to which was added 50  $\mu$ L β-mercaptoethanol 0.2% at final concentration, 0.1 mg (3.2 U) of proteinase K (modified method of Quinta et al., 2004), and incubated at 65°C for 12 h. The solution was then centrifuged at 12,000 g (10 min), the supernatant was extracted with chloroform–isoamylic alcohol and precipitated with ethanol. The pellets were finally re-suspended in TE buffer.

### Method using TRITON/CI chloroform–isoamyl alcohol (EX6)

The muscle tissue (0.1 g) was homogenized in 500  $\mu$ L of pH 8.0-buffer B (2 mM EDTA,10 mM Tris–HCl), to which was added 50  $\mu$ L β-mercaptoethanol 0.2% at final concentration and 50  $\mu$ L Triton X-100 at 1% final concentration and 0.1 mg (3.2 U) of proteinase K (Besbes et al., 2011). The mixture was incubated at 65°C for 12 h, and then centrifuged at at 12,000 g (4°C-10 min). The suspended DNA was mixed vigorously with an equal volume of chloroform–isoamylic alcohol and then centrifuged at 12,000g (10 min). For DNA precipitation, ethanol and sodium acetate were added to the supernatant, the pellets were then resuspended in 50  $\mu$ L TE buffer.

#### Method using DNeasy tissue kit (EX7)

DNA was extracted by means of a commercial kit (DNeasy Tissue kit, Qiagen, Darmstadt, Germany), as previously described (Prado et al., 2005), 25 mg of sample were homogenized with 180  $\mu$ l of lysis buffer ATL and 20  $\mu$ l proteinase K. The homogenate was incubated during 3 h at 56°C and then Vortex for 15 s. 200  $\mu$ l Buffer AL were then added to the sample, and thoroughly mixed by vortexing. The supernatant, obtained after centrifugation at 6000 g for 5 min at room temperature, was added to 200  $\mu$ l ethanol. The mixture (400  $\mu$ l) was then eluated on the DNeasy Mini spin column where DNA was bonded to the silica-gel membrane during the centrifugation steps. The bounded DNA-DNeasy Mini spin membrane was washed using two different washing buffers, and then eluted using the elution buffer.

#### **Evaluation of DNA quantification and purity**

The quantification of DNA was measured using a spectrophotometer (absorbance at 260 nm) and its purity was checked using the ratio A260/A280 as described by Sambrock et al. (1989).

#### **PCR** Amplification

The purified DNA was used for PCR amplification. As processing cause DNA degradation, the primers UnivPF (50-GACAATACCAATTCGTAAATC) and UnivPR (50 -CCATAATAAATACCTCGTCC-30) designed using Mega software and based on the highly conserved sequence among the target species and were used for this study, allowed the amplification of a 315 bp fragment. A final reactional volume of 5 µL was used for PCR amplification. For each reaction, 100 ng of extracted DNA template were mixed with 5 µL of 5xPCR buffer, 2 µL of dNTP mix (10 mM of each); 50 pmol of each primer, 3 mM of MgCl2 and 5 U of Taq DNA polymerase (Go-Taq, Promega, France). The PCR conditions used for amplification passed by a step of denaturation (94°C for 5 min), ensued by 35 cycles consisting of 94°C for 30s, 55°C for 30s, 72°C for 40s, and finally by an extension step (72°C for 7 min). Eelectrophoresis (2% agarose gel, TBE buffer containing 1 mM EDTA, 45 mM Tris borate, pH 8, and 0.2  $\mu$ g/mL ethidium bromide) was used to check DNA amplification. Before sequencing, the PCR products were visualised via ultraviolet trans-illumination. A commercial 100pb DNA ladder (Fermentas, France) was used to compare the size of the resulting DNA fragments.

#### **Statistical Analysis**

All results were presented as mean  $\pm$  standard deviation (n=6) and statistical analyses were performed using ANOVA with an SPSS Statistical Software System 11.01 (SPSS Chicago, IL, USA). Eventual significant differences among mean values (5% level) were checked using Tukey's test. Data were treated by partial least squares (PLS1) regression using multivariate statistical software (Unscrambler version 9.8, CAMO Software AS, Oslo, Norway), and were subjected to Leverage correction and cross validation. In order to compensate for the different scales of the variables, they were weighted with the inverse of the standard deviation of all the objects. The influence of the studied variables on PCR amplification was examined using Martens Uncertainty test.

#### **Results and Discussion**

## Evaluation of DNA yield, purity and amplificability

The different samples used to evaluate the DNA yield and purity, included the tissue of fresh shrimp and that of moderately (frozen) and highly processed (breaded and cooked) products. Gel electrophoresis and spectrophotometry analysis revealed significant differences (p<0.05) in DNA yield and purity according to the procedure of extraction and sample matrix (Table1). In the process of DNA extraction, 3 steps are necessary including breaking cell walls, nucleic material release from the nuclei into the media and DNA precipitation. In this study, the highest precipitated amount of DNA was obtained following the PCI/SDS protocol and that, independently of fresh or processed shrimp (Table 1). However, adding NaCl to the PCI/SDS medium induced a significant decrease in DNA levels in all analysed sample. Such result could be due to the inadequate concentration used in this experiment as NaCl has either dispersive or binding effects on DNA fragment causing its precipitation along with proteins (Bauer et al., 2013, Besbes et al., 2011). However, no optimisation of NaCl concentration was performed in this study as sufficient quantity of DNA level with good purity (A260/A280 ratio >1.81) was obtained in all cases.

Significant lower DNA contents were obtained for cooked and breaded shrimp when using specific method such as the Chelex or DNeasy Tissue Kit procedure **(Table 1)**. Moreover for such processed samples, both methods yielded poor DNA quality (A260/A280 <1.21) which may be attributed to interferences with substances such as protein and flavouring aromatic compounds added to shrimp preparation.

#### Journal abbreviation: J FisheriesSciences.com

 Table 1: Comparison of DNA yield, purity, PCR positive amplification for several foodstuffs of shrimp extracted with six different methods studied.

Kit Fast DNA MP         32.9 ± 1.47         2.08.00125         +           Penaeus keraihurus         IOTISDS MCI         277.41 ± 9.41%         1.97 ± 0.02         +           PCINDS/NACI         277.41 ± 9.41%         1.97 ± 0.02         +           CHELEX         308.1 ± 16.66%         1.87 ± 0.06         +           CTAB         617.5 ± 39.56%         1.87 ± 0.06         +           TRITON         635.8 ± 37.27         1.86 ± 0.01         +           Parapenaeus Jongirostris         Q2.9 ± 1.82         1.88 ± 0.01         +           PCISDS MARCI         228.26 ± 8.17%         1.97 ± 0.02         +           PCISDSNACI         228.26 ± 8.17%         1.97 ± 0.02         +           C1AB         40.01 ± 21.93*         2.01 ± 0.015         +           C1AB         40.01 ± 21.93*         2.01 ± 0.015         +           TRITON         601.75 ± 28.49*         1.88 ± 0.005         +           TRITON         617.5 ± 28.49*         1.98 ± 0.015         +           PCISDSNACI         224.4         1.82 ± 0.015         +           TRITON         625.35 ± 30.05*         1.91 ± 0.045         +           C1AB         558.4 ± 2.069         1.99 ± 0.005         + <t< th=""><th>Sample</th><th><b>Extraction methods</b></th><th>DNA (ng/µl extract)</th><th><math>A_{260}/A_{280}</math></th><th>PCR</th></t<>	Sample	<b>Extraction methods</b>	DNA (ng/µl extract)	$A_{260}/A_{280}$	PCR
Peraeus kerathurus         PCISDS         1601.78 ± 15.48 <sup>b</sup> 2.05 ± 0.03         + + PCISDSNaCI           Peraeus kerathurus         CHELEX         308.1 ± 16.66 <sup>c</sup> 1.87 ± 0.06         + - CTAB         617.5 ± 39.56 <sup>c</sup> 1.87 ± 0.06         + - DReasy Tissue kit         92.9 ± 1.52 <sup>t</sup> 1.88 ± 0.01         + - - -           FRIF         TRITON         635.8 ± 37.2 <sup>c</sup> 1.88 ± 0.01         + - -         -         + -         -         -         -         -         -         -         -         + -         -	-	Kit Fast DNA MP	$32.9 \pm 1.47^{a}$	$2.08 \pm 0.0125$	+
Penaeus kerathurus         PCI/SDS/NaCl         277.41 ± 9.41*         1.97 ± 0.02         + +           CHRLEX         308.1 ± 16.66*         1.87 ± 0.06         +           CHRLEX         308.1 ± 16.66*         1.87 ± 0.06         +           TRUION         635.8 ± 37.2*         1.86 ± 0.01         +           Dreasy Tissue kit         92.9 ± 1.52*         1.88 ± 0.01         +           Parapenaeus longirostris         Kit Fast DNA MP         41.6 ± 1.8 *         2.08 ± 0.062         +           PCISDS/NaCl         228.26 ± 8.17*         1.97 ± 0.02         +         +           PCISDS/NaCl         228.26 ± 8.17*         1.97 ± 0.02         +         +           CTAB         40.01 ± 2.1.93*         2.01 ± 0.015         +         +           TRITON         601.75 ± 28.49*         1.98 ± 0.005         +         +           DNexsy Tissue kit         1411.40 ± 13.10*         2.01 ± 0.015         +           PCI/SDS/NaCl         284.5 ± 17.03*         2.05 ± 0.045         +           CHELEX         22.03 ± 3.01*         1.91 ± 0.015         +           PCI/SDS/NaCl         284.5 ± 17.03*         2.05 ± 0.045         +           CHELEX         22.03 ± 3.01*         1.91 ± 0.015 <t< td=""><td>PCI/SDS</td><td><math>1691.78 \pm 15.48^{\rm b}</math></td><td><math display="block">2.05\pm0.03</math></td><td>+</td></t<>		PCI/SDS	$1691.78 \pm 15.48^{\rm b}$	$2.05\pm0.03$	+
Penaeus kerathurus         CHELEX         308.1 ± 16.66°         1.87 ± 0.06         +           CTAB         617.5 ± 39.56°         1.87 ± 0.06         +           TRITON         635.8 ± 37.2°         1.86 ± 0.01         +           DNessy Tissue kit         92.9 ± 1.52′         1.88 ± 0.01         +           KIF fast DNA MP         41.6 ± 1.8°         2.08 ± 0.062         +           PCISDSNaCL         128.26 ± 8.17°         1.97 ± 0.02         +           PCISDSNacL         228.26 ± 8.17°         1.97 ± 0.02         +           CTAB         440.1 ± 21.93 <sup>4</sup> 2.01 ± 0.015         +           TRITON         60.175 ± 28.49°         1.98 ± 0.005         +           TRITON         60.175 ± 28.49°         1.98 ± 0.005         +           TRITON         60.175 ± 28.49°         1.98 ± 0.005         +           TRITON         60.175 ± 28.49°         1.98 ± 0.015         +           PCISDSNaCL         284.5 ± 17.03°         2.05 ± 0.045         +           TRITON         625.35 ± 30.05°         1.98 ± 0.01         +           PCISDSNaCL         284.5 ± 17.03°         2.05 ± 0.045         +           TRITON         625.35 ± 30.05°         1.98 ± 0.01         + <td></td> <td>PCI/SDS/NaCl</td> <td>277.41 ± 9.41°</td> <td><math>1.97 \pm 0.02</math></td> <td>+</td>		PCI/SDS/NaCl	277.41 ± 9.41°	$1.97 \pm 0.02$	+
CTAB         617.5 ± 39.50*         1.87 ± 0.065         +           TRITON         635.8 ± 37.2*         1.86 ± 0.01         +           DNeasy Tissue kit         92.9 ± 1.52*         1.88 ± 0.01         +           Kit Fast DNA MP         41.6 ± 1.8*         2.08 ± 0.062         +           PCUSDS         1486.46 ± 14.4%         2.01 ± 0.035         +           PCUSDSNaCI         228.26 ± 8.17*         1.97 ± 0.02         +           CHELEX         290.1 ± 32.4*         1.82 ± 0.015         +           CTAB         44.0 1 ± 2.1.93*         2.01 ± 0.015         +           CTAB         44.0 1 ± 2.1.93*         2.01 ± 0.015         +           CTAB         44.0 1 ± 0.1.93*         2.01 ± 0.015         +           DNeasy Tissue kit         84.9 ± 2.7*         1.72 ± 0.02         +           Kit Fast DNA         41.1 ± 0.50*         1.86 ± 0.015         +           PCISDS         1411.46 ± 13.10*         2.01 ± 0.015         +           PCISDSNCI         228.45 ± 17.03*         2.03 ± 0.045         -           CTAB         58.8 ± 22.69*         1.98 ± 0.01         +           TRITON         652.3 ± 10.05         +            PCISDSNACI <t< td=""><td>Penaeus kerathurus</td><td>CHELEX</td><td><math>308.1 \pm 16.66^{\circ}</math></td><td><math display="block">1.87\pm0.06</math></td><td>+</td></t<>	Penaeus kerathurus	CHELEX	$308.1 \pm 16.66^{\circ}$	$1.87\pm0.06$	+
TRITON         653 k ± 37 2*         1.86 ± 0.01         +           DNeasy Tissue kit         92.9 ± 1.52 f         1.88 ± 0.01         +           Kit Fast DNA MP         41.6 ± 1.83         2.08 ± 0.062         +           PCISDSN         1486.46 ± 1.4.44*         2.01 ± 0.035         +           PCISDSNC1         22.82 d ± 8.17*         1.97 ± 0.02         +           PCISDSNC1         22.82 d ± 8.17*         1.97 ± 0.02         +           CHELEX         290.1 ± 32.4*         1.82 ± 0.015         +           CHBA         440.1 ± 2.193*         2.01 ± 0.015         +           TRITON         601.75 ± 28.49*         1.98 ± 0.005         +           TRISON         611.1 ± 0.50*         1.86 ± 0.015         +           Kit Fast DNA         41.1 ± 0.50*         1.86 ± 0.015         +           PCI/SDS         1411.46 ± 13.10*         2.01 ± 0.015         +           PCI/SDS/NACI         28.45 ± 1.7.03*         2.05 ± 0.045         +           TRITON         625.35 ± 30.05*         1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4*         1.88 ± 0.04         +           Kit Fast DNA MP         41.7 ± 0.71*         1.87 ± 0.01         +		СТАВ	$617.5 \pm 39.56^{\rm e}$	$1.87 \pm 0.065$	+
DNeasy Tissue kit         92.9 ± 1.52*         1.88 ± 0.01         +           Kit Fast DNA MP         41.6 ± 1.8*         2.08 ± 0.062         +           Parapenaeus longirostris         I486.46 ± 14.44 <sup>3</sup> 2.01 ± 0.035         +           PCI/SDS         1486.46 ± 14.44 <sup>3</sup> 2.01 ± 0.035         +           PCI/SDS/NaCl         228.26 ± 8.17*         1.97 ± 0.02         +           CIAB         440.1 ± 21.93*         2.01 ± 0.015         +           CIAB         440.1 ± 21.93*         2.01 ± 0.015         +           TRITON         601.75 ± 28.49*         1.98 ± 0.005         +           DNeasy Tissue kit         84.9 ± 2.7*1         1.72 ± 0.02         +           Kit Fast DNA         41.1 ± 0.50*         1.88 ± 0.015         +           PCUSDS         1411.46 ± 13.10*         2.01 ± 0.015         +           PCUSDS/NaCl         284.5 ± 17.03*         2.05 ± 0.045         -           CIAB         585.8 ± 22.60*         1.99 ± 0.015         +           DNeasy Tissue kit         105.6 7 ± 1.4*         1.88 ± 0.01         +           DNeasy Tissue kit         107.9 ± 2.65.2*         1.98 ± 0.015         +           PCUSDS/NaCl         233.2 ± 17.3*         1.99 ± 0.015		TRITON	$635.8 \pm 37.2^{\circ}$	$1.86 \pm 0.01$	+
Kit Fast DNA MP         41.6 ± 1.8*         2.08 ± 0.062         +           Parapenaeus longirostris         PCUSDS         1486.46 ± 14.44*         2.01 ± 0.035         +           PCUSDS/NACI         228.26 ± 8.17*         1.97 ± 0.02         +           CHELEX         290.1 ± 32.4*         1.82 ± 0.015         +           CTAB         440.1 ± 21.93 <sup>3</sup> 2.01 ± 0.015         +           TRITON         6601.75 ± 28.49*         1.98 ± 0.005         +           DNeasy Tissue kit         84.9 ± 2.7 f         1.72 ± 0.02         +           Kit Fast DNA         41.1 ± 0.50*         1.86 ± 0.015         -           PCUSDS         1411.46 ± 13.10*         2.01 ± 0.015         +           PCUSDS         1411.46 ± 13.10*         2.01 ± 0.015         +           CTAB         583.8 ± 2.69 <sup>3</sup> 1.98 ± 0.01         +           CTAB         583.8 ± 2.69 <sup>3</sup> 1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4*         1.88 ± 0.04         +           CTAB         583.8 ± 2.69 <sup>3</sup> 1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4*         1.88 ± 0.04         +           CTAB         41.23 ± 13.57*         1.99 ± 0.005		DNeasy Tissue kit	$92.9 \pm 1.52$ f	$1.88 \pm 0.01$	+
PCUSDS         1486.46 ± 14.44°         2.01 ± 0.035         +           PCUSDS/NaC1         228.26 ± 8.17°         1.97 ± 0.02         +           CHELEX         290.1 ± 32.4°         1.82 ± 0.015         +           CTAB         440.1 ± 21.93 <sup>4</sup> 2.01 ± 0.015         +           TRITON         601.75 ± 28.49°         1.98 ± 0.005         +           DNeasy Tissue kit         84.9 ± 2.7 °         1.72 ± 0.02         +           Kit Fast DNA         41.1 ± 0.50°         1.88 ± 0.015         -           PCUSDS         1411.46 ± 13.10°         2.01 ± 0.015         +           PCUSDS/NaC1         284.5 ± 17.03°         2.05 ± 0.045         -           CTAB         585.8 ± 22.69 <sup>4</sup> 1.90 ± 0.055         +           TRITON         625.35 ± 30.05 <sup>4</sup> 1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4°         1.88 ± 0.04         +           TRITON         625.35 ± 26.52 <sup>h</sup> 1.98 ± 0.015         +           PCI/SDS/NaC1         233.2 ± 17.34 <sup>k</sup> 195 ± 0.015         +           PCI/SDS/NaC1         233.2 ± 17.34 <sup>k</sup> 195 ± 0.015         +           PCI/SDS/NaC1         233.2 ± 17.34 <sup>k</sup> 195 ± 0.015         +		Kit Fast DNA MP	$41.6 \pm 1.8^{a}$	$2.08 \pm 0.062$	+
Parapenaeus longirostris         PCL'SDS/NaC1         228 26 ± 8.17*         1.97 ± 0.02         +           CHE         290.1 ± 32.4*         1.82 ± 0.015         +           CTAB         440.1 ± 2.193*         2.01 ± 0.015         +           TRITON         601.75 ± 28.49*         1.98 ± 0.005         +           DNeasy Tissue kit         84.9 ± 2.7*         1.72 ± 0.02         +           Kit Fast DNA         41.1 ± 0.50*         1.86 ± 0.015         -           PCUSDSNaC1         228.4 ± 17.03*         2.05 ± 0.045         -           PCUSDSNaC1         228.4 ± 17.03*         2.05 ± 0.045         -           PCUSDSNaC1         228.4 ± 17.03*         2.05 ± 0.045         +           PCUSDSNaC1         228.4 ± 17.03*         2.05 ± 0.045         +           CHELEX         220.3 ± 0.05*         +         +           CHELEX         109.0 ± 0.055         +         +           TRITON         62.53 ± 3.005*         1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4*         1.88 ± 0.04         +           Kit Fast DNA MP         41.7 ± 0.71*         1.87 ± 0.06         -           COLVISDS/NaC1         233.2 ± 17.34*         1.95 ± 0.015         +		PCI/SDS	$1486.46 \pm 14.44^{b}$	$2.01 \pm 0.035$	+
Parapenaeus longirostris         CHELEX         290.1 ± 32.4*         1.82 ± 0.015         +           CTAB         440.1 ± 21.93 <sup>d</sup> 2.01 ± 0.015         +           DNeasy Tissue kit         84.9 ± 2.7 f         1.78 ± 0.005         +           Metapenaeus monoceros         CTXBS         411.14 ± 0.50 <sup>d</sup> 1.86 ± 0.015         -           PC/SDS         1411.46 ± 13.10 <sup>b</sup> 2.01 ± 0.015         +         -           PC/SDS/NaC1         284.5 ± 17.03 <sup>s</sup> 2.05 ± 0.045         -         -           CHELEX         220.3 ± 3.61 <sup>s</sup> 1.91 ± 0.045         +         -         -           PC/SDS/NaC1         284.5 ± 17.03 <sup>s</sup> 2.05 ± 0.045         +         -         -           DNeasy Tissue kit         105.67 ± 1.4 <sup>s</sup> 1.88 ± 0.04         +         +         -	Parapenaeus longirostris	PCI/SDS/NaCl	$228.26 \pm 8.17^{\circ}$	$1.97 \pm 0.02$	+
$Frozen peeled precooked shrimp \\ Frozen peeled Shrimp \\ Frozen Peele Shrimp \\ Frozen Peele Shrimp \\ Frozen Peele Shrimp \\ Frozen Peeleele Shrimp \\ Frozen Peele Shrimp \\ Frozen Peele Shrim$		CHELEX	290.1 ± 32.4°	$1.82 \pm 0.015$	+
TRITON         601.75 ± 28.49*         1.98 ± 0.005         +           DNeasy Tissue kit         84.9 ± 2.7 '         1.72 ± 0.02         +           Kit Fast DNA         41.1 ± 0.50*         1.86 ± 0.015         -           PCI/SDS         1411,46 ± 13,10*         2.01 ± 0.015         +           PCI/SDS/NaCl         284.5 ± 1.7.03*         2.05 ± 0.045         -           CHELEX         22.0.3 ± 3.61*         1.91 ± 0.045         +           CTAB         585.8 ± 22.694         1.90 ± 0.055         +           TRUTON         625.35 ± 30.054         1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4*         1.88 ± 0.04         +           Kit Fast DNA MP         41.7 ± 0.71*         1.87 ± 0.06         -           PCI/SDS/NaCl         23.5 ± 17.34*         1.95 ± 0.015         +           PCI/SDS/NaCl         23.5 ± 17.34*         1.95 ± 0.015         +           CTAB         412.3 ± 13.57*         1.99 ± 0.005         +           TRITON         58.64 ± 9.98*         1.84 ± 0.02         +           DNeasy Tissue kit         84.17 ± 0.6 4         1.82 ± 0.01         +           Kit Fast DNA MP         12.5 ± 7.24*         1.61 ± 0.01         -		СТАВ	$440.1 \pm 21.93^{d}$	$2.01 \pm 0.015$	+
DNeasy Tissue kit         84.9 ± 2.7 f         1.72 ± 0.02         +           Kit Fast DNA         41.1 ± 0.50°         1.86 ± 0.015         -           PCI/SDS         1411,46 ± 13,10°         2.01 ± 0.015         +           PCI/SDS         1244,5 ± 17.03°         2.05 ± 0.045         -           CHELEX         220.3 ± 3.61°         1.91 ± 0.045         +           CTAB         585,8 ± 22.69 <sup>4</sup> 1.90 ± 0.055         +           TRITON         625.35 ± 30.05 <sup>4</sup> 1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4°         1.88 ± 0.04         +           Kit Fast DNA MP         41.7 ± 0.71°         1.87 ± 0.06         -           PCI/SDS         1679.5 ± 26.52 <sup>b</sup> 1.98 ± 0.015         +           PCI/SDS         1679.5 ± 26.52 <sup>b</sup> 1.98 ± 0.015         +           PCI/SDSNACI         23.3.2 ± 17.34°         1.95 ± 0.05         +           TRITON         588.61 ± 29.98°         1.84 ± 0.02         -           CTAB         412.3 ± 13.57°         1.99 ± 0.005         +           TRITON         588.61 ± 29.98°         1.61 ± 0.01         -           PCI/SDS         1015.6 ± 4.95 <sup>6</sup> 2.02 ± 0.015         -		TRITON	$601.75 \pm 28.49^{\circ}$	$1.98 \pm 0.005$	+
Kit Fast DNA         41.1 ± 0.50°         1.86 ± 0.015         -           PCV/SDS         1411,46 ± 13,10°         2.01 ± 0.015         +           PCVSDS/NaCI         224.5 ± 17.03°         2.05 ± 0.045         -           CHELEX         220.3 ± 3.61°         1.91 ± 0.045         +           CTAB         585.8 ± 22.69°         1.90 ± 0.055         +           TRITON         662.35 ± 30.05°         1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4°         1.88 ± 0.04         +           Kit Fast DNA MP         41.7 ± 0.71°         1.88 ± 0.04         +           PC/SDS/NaCI         233.2 ± 17.34°         1.95 ± 0.015         +           PC/SDS/NaCI         233.2 ± 17.34°         1.95 ± 0.015         +           PC/SDS/NaCI         233.2 ± 17.34°         1.95 ± 0.015         +           PC/SDS/NaCI         233.2 ± 17.34°         1.99 ± 0.002         -           CTAB         412.3 ± 13.57°         1.99 ± 0.005         +           TRITON         58.86 ± 29.98°         1.84 ± 0.025         +           DNeasy Tissue kit         84.17 ± 0.64         1.82 ± 0.01         +           Reito precooked shrimp         CHELEX         53.03 ± 11.64         1.61 ± 0.01		DNeasy Tissue kit	$84.9 \pm 2.7$ f	$1.72 \pm 0.02$	+
Metapenaeus monoceros         PCUSDS         1411,46 ± 13,10°         2.01 ± 0.015         +           Metapenaeus monoceros         CHELEX         220.3 ± 0.10°         2.05 ± 0.045         -           CHELEX         220.3 ± 3.61°         1.91 ± 0.045         +         +           CTAB         585.8 ± 22.69 <sup>4</sup> 1.90 ± 0.055         +           TRITON         625.35 ± 30.05°         1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4°         1.88 ± 0.04         +           Kit Fast DNA MP         41.7 ± 0.71°         1.87 ± 0.06         -           PCI/SDS         1679.5 ± 26.52°         1.98 ± 0.015         +           PCI/SDSNaC1         233.2 ± 17.34°         1.95 ± 0.015         +           PCI/SDSNaC1         233.2 ± 17.34°         1.95 ± 0.015         +           PCI/SDSNaC1         233.2 ± 17.34°         1.95 ± 0.015         +           TRITON         588.61 ± 29.98°         1.84 ± 0.025         +           TRITON         588.61 ± 29.98°         1.84 ± 0.025         +           DNeasy Tissue kit         84.17 ± 0.64°         1.82 ± 0.01         +           PCI/SDS         1015.6 ± 4.95°         2.01 ± 0.005         +           PCI/SDS/NaC1		Kit Fast DNA	$41.1 \pm 0.50^{a}$	$1.86 \pm 0.015$	-
Metapenaeus monoceros         PCL/SDS/NaCl         284.5 ± 17.03°         2.05 ± 0.045         -           CHEL.EX         220.3 ± 3.61°         1.91 ± 0.045         +           CTAB         855.8 ± 22.694         1.90 ± 0.055         +           TRITON         625.35 ± 30.054         1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4°         1.88 ± 0.04         +           Kit Fast DNA MIP         41.7 ± 0.71*         1.87 ± 0.06         -           PCI/SDS         1679.5 ± 26.52 <sup>b</sup> 1.98 ± 0.015         +           PCL/SDS/NaCI         233.2 ± 17.34*         1.95 ± 0.015         +           PCL/SDS/NaCI         233.2 ± 17.34*         1.95 ± 0.015         +           CTAB         412.3 ± 13.57°         1.99 ± 0.005         +           TRITON         588.61 ± 29.98'         1.84 ± 0.025         +           DNeasy Tissue kit         84.17 ± 0.61         6.182 ± 0.01         +           PCI/SDS         1015.6 ± 4.95 <sup>b</sup> 2.01 ± 0.005         +           PCI/SDS/NaCI         130.86 ± 7.58 <sup>c</sup> 2.02 ± 0.015         -           Coked shrimp         CHELEX         53.03 ± 11.16 <sup>4</sup> 1.13 ± 0.02         -           CHELEX         53.03 ± 11.50 <sup></sup>	-	PCI/SDS	$1411,46 \pm 13,10^{b}$	$2.01 \pm 0.015$	+
Metapenaeus monoceros         CHELEX         220.3 ± 3.61 <sup>e</sup> 1.91 ± 0.045         +           CTAB         585.8 ± 22.69 <sup>d</sup> 1.90 ± 0.055         +           TRITON         625.35 ± 30.05 <sup>d</sup> 1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4 <sup>e</sup> 1.88 ± 0.04         +           Kit Fast DNA MP         41.7 ± 0.71 <sup>s</sup> 1.87 ± 0.06         -           PCI/SDS         1679.5 ± 26.52 <sup>b</sup> 1.98 ± 0.015         +           PCI/SDS/NaCI         233.2 ± 17.34 <sup>e</sup> 1.95 ± 0.015         +           CTAB         412.3 ± 13.57 <sup>c</sup> 1.99 ± 0.02         -           CTAB         412.3 ± 13.57 <sup>c</sup> 1.99 ± 0.005         +           TRITON         588.61 ± 29.98 <sup>c</sup> 1.84 ± 0.025         +           DNeasy Tissue kit         84.17 ± 0.6 <sup>4</sup> 1.82 ± 0.01         +           Kit Fast DNA MP         12.5 ± 7.24 <sup>s</sup> 1.61 ± 0.01         -           PCI/SDSNaCI         13.08 ± 7.58 <sup>s</sup> 2.01 ± 0.005         +           PCI/SDSNaCI         13.08 ± 7.58 <sup>s</sup> 2.00 ± 0.02         +           TRITON         328.1 ± 11.96 <sup>c</sup> 2.00 ± 0.02         +           Cooked shrimp         CTAB         450.2		PCI/SDS/NaCl	$284.5 \pm 17.03^{\circ}$	$2.05 \pm 0.045$	-
	Metapenaeus monoceros	CHELEX	$220.3 \pm 3.61^{\circ}$	$1.91 \pm 0.045$	+
TRITON         625.35 ± 30.05 <sup>4</sup> 1.98 ± 0.01         +           DNeasy Tissue kit         105.67 ± 1.4°         1.88 ± 0.04         +           Kit Fast DNA MP         41.7 ± 0.71 <sup>a</sup> 1.87 ± 0.06         -           PCI/SDS         1679.5 ± 26.52°         1.98 ± 0.015         +           PCI/SDS/NaCI         233.2 ± 17.34°         1.95 ± 0.015         +           PCI/SDS/NaCI         233.2 ± 17.34°         1.99 ± 0.02         -           CTAB         412.3 ± 13.57°         1.99 ± 0.005         +           TRITON         588.61 ± 29.98°         1.84 ± 0.023         +           DNeasy Tissue kit         84.17 ± 0.6 <sup>d</sup> 1.82 ± 0.01         +           Kit Fast DNA MP         12.5 ± 7.24°         1.61 ± 0.01         -           PC/SDS/NaCI         130.86 ± 7,58°         2.02 ± 0.005         +           PC/SDS/NaCI         130.86 ± 7,58°         2.02 ± 0.015         -           CTAB         450.26 ± 1.246°         2.01 ± 0.02         -           TRITON         328.1 ± 11.96°         2.00 ± 0.02         +           DNeasy Tissue kit         87.17 ± 3 ľ         1.71 ± 0.01         +           Kit Fast DNA MP         10.8 ± 0.16°         1.69 ± 0.01         -		CTAB	$585.8 \pm 22.69^{d}$	$1.90 \pm 0.055$	+
DNeasy Tissue kit $105.67 \pm 1.4^{\circ}$ $1.88 \pm 0.04$ $+$ Kit Fast DNA MP $41.7 \pm 0.71^{\circ}$ $1.87 \pm 0.06$ $-$ PCI/SDS $1679.5 \pm 60.26^{\circ}$ $1.98 \pm 0.015$ $+$ PC/SDS/NaCI $233.2 \pm 17.34^{\circ}$ $1.95 \pm 0.015$ $+$ CHELEX $89.0 \pm 2.67^{\circ}$ $1.39 \pm 0.02$ $-$ CTAB $412.3 \pm 13.57^{\circ}$ $1.99 \pm 0.005$ $+$ TRITON $588.61 \pm 29.98'$ $1.84 \pm 0.025$ $+$ DNeasy Tissue kit $84.17 \pm 0.6^{\circ}$ $1.82 \pm 0.01$ $+$ Kit Fast DNA MP $12.5 \pm 7.24^{\circ}$ $1.61 \pm 0.01$ $-$ PCU/SDS/NaCI $130.86 \pm 7.58^{\circ}$ $2.02 \pm 0.015$ $-$ Cooked shrimpCHELEX $53.03 \pm 11.16^{\circ}$ $1.13 \pm 0.02$ $-$ Cooked shrimpCHELEX $53.03 \pm 11.16^{\circ}$ $2.00 \pm 0.02$ $+$ PCU/SDS/NaCI $108.20.66 \pm 12.46^{\circ}$ $2.00 \pm 0.02$ $+$ TRITON $328.1 \pm 11.96^{\circ}$ $2.00 \pm 0.02$ $+$ DNeasy Tissue kit $87.17 \pm 3^{\circ}$ $1.71 \pm 0.01$ $+$ Kit Fast DNA MP $10.8 \pm 0.16^{\circ}$ $1.69 \pm 0.01$ $-$ CTAB $450.26 \pm 12.46^{\circ}$ $2.00 \pm 0.02$ $+$ DNeasy Tissue kit $87.17 \pm 3^{\circ}$ $1.71 \pm 0.01$ $+$ Kit Fast DNA MP $10.8 \pm 0.16^{\circ}$ $1.69 \pm 0.01$ $-$ CCAB $62.3 \pm 12.0^{\circ}$ $1.97 \pm 0.03$ $+$ TRITON $391.9 \pm 17.8^{\circ}$ $1.73 \pm 0.01$ $-$ CDNeasy Tissue kit $64.13 \pm 8^{\circ}$ $1.41 \pm 0.01$ $+$ Kit Fast		TRITON	$625.35 \pm 30.05^{d}$	$1.98 \pm 0.01$	+
Kit Fast DNA MP $41.7 \pm 0.71^{3}$ $1.87 \pm 0.06$ -PCl/SDS $1679.5 \pm 26.52^{5}$ $1.98 \pm 0.015$ +PCl/SDS/NaCl $233.2 \pm 17.34^{c}$ $1.95 \pm 0.015$ +PCl/SDS/NaCl $233.2 \pm 17.34^{c}$ $1.95 \pm 0.015$ +CTAB $412.3 \pm 13.57^{a}$ $1.99 \pm 0.005$ +TRITON $588.61 \pm 29.98^{t}$ $1.84 \pm 0.025$ +DNeasy Tissue kit $84.17 \pm 0.6^{-d}$ $1.82 \pm 0.01$ +Kit Fast DNA MP $12.5 \pm 7.24^{a}$ $1.61 \pm 0.01$ -PCl/SDS $1015.6 \pm 4.95^{b}$ $2.01 \pm 0.005$ +PCl/SDS/NaCl $130.86 \pm 7.58^{-c}$ $2.02 \pm 0.015$ -Cooked shrimpCHELEX $53.03 \pm 11.16^{d}$ $1.13 \pm 0.02$ -CTAB $450.26 \pm 12.46^{c}$ $2.01 \pm 0.02$ -TRITON $328.1 \pm 11.96^{c}$ $2.00 \pm 0.02$ +DNeasy Tissue kit $87.17 \pm 3^{-c}$ $1.71 \pm 0.01$ +Kit Fast DNA MP $10.8 \pm 0.16^{a}$ $1.69 \pm 0.01$ -PCL/SDS/NaCl $59.51 \pm 7.39^{c}$ $1.73 \pm 0.01$ -CTAB $62.3 \pm 1.20^{c}$ $1.97 \pm 0.03$ +TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ +DNeasy Tissue kit $64.13 \pm 8^{-c}$ $1.41 \pm 0.01$ +Kit Fast DNA MP $19.4 \pm 0.3^{3} \pm 1.71 \pm 0.03$ +TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ +DNeasy Tissue kit $64.13 \pm 8^{-c}$ $1.41 \pm 0.01$ +CTAB $62.3 \pm 1.20^{c}$ $1.97 \pm 0.03$ + <t< td=""><td>DNeasy Tissue kit</td><td>105.67 ± 1.4 °</td><td><math>1.88 \pm 0.04</math></td><td>+</td></t<>		DNeasy Tissue kit	105.67 ± 1.4 °	$1.88 \pm 0.04$	+
Plesionika heterocarpusPCL/SDS $1679.5 \pm 26.52^{b}$ $1.98 \pm 0.015$ $+$ PCl/SDS/NaCl $233.2 \pm 17.34^{a}$ $1.95 \pm 0.015$ $+$ CHELEX $89.0 \pm 2.67^{d}$ $1.39 \pm 0.02$ $-$ CTAB $412.3 \pm 13.57^{e}$ $1.99 \pm 0.005$ $+$ TRITON $588.61 \pm 29.98^{t}$ $1.84 \pm 0.025$ $+$ DNeasy Tissue kit $84.17 \pm 0.6^{4}$ $1.82 \pm 0.01$ $+$ Kit Fast DNA MP $12.5 \pm 7.24^{s}$ $1.61 \pm 0.01$ $-$ PCL/SDS $1015.6 \pm 4.95^{b}$ $2.01 \pm 0.005$ $+$ PCU/SDS/NaCl $130.86 \pm 7.58^{e}$ $2.02 \pm 0.015$ $-$ Cooked shrimpCHELEX $53.03 \pm 11.16^{d}$ $1.13 \pm 0.02$ $-$ CTAB $450.26 \pm 12.46^{e}$ $2.01 \pm 0.00$ $-$ TRITON $328.1 \pm 11.96^{e}$ $2.00 \pm 0.02$ $+$ DNeasy Tissue kit $87.17 \pm 3^{t}$ $1.71 \pm 0.01$ $+$ Frozen peeled precooked shrimpCHELEX $53.03 \pm 11.16^{d}$ $1.82 \pm 0$ $+$ PCL/SDS $1090.03 \pm 11.05^{b}$ $1.82 \pm 0$ $+$ $-$ PCL/SDS/NaCl $59.51 \pm 7.39^{e}$ $1.73 \pm 0.01$ $-$ CTAB $62.3 \pm 1.20^{e}$ $1.97 \pm 0.03$ $+$ TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ $+$ DNeasy Tissue kit $64.13 \pm 8^{c}$ $1.41 \pm 0.01$ $+$ CYSDS/NaCl $47.35 \pm 8.70^{c}$ $1.36 \pm 0.005$ $-$ CTAB $62.3 \pm 1.20^{e}$ $1.32 \pm 0.01$ $+$ TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ $+$	Plesionika heterocarpus	Kit Fast DNA MP	$41.7 \pm 0.71^{a}$	$1.87 \pm 0.06$	-
Plesionika heterocarpus $PCU/SDS/NaCl$ $233.2 \pm 17.34^{\circ}$ $1.95 \pm 0.015$ $+$ CHELEX $89.0 \pm 2.67^{d}$ $1.39 \pm 0.02$ $-$ CTAB $412.3 \pm 13.57^{\circ}$ $1.99 \pm 0.005$ $+$ TRITON $588.61 \pm 29.98^{t}$ $1.84 \pm 0.025$ $+$ DNeasy Tissue kit $84.17 \pm 0.6^{d}$ $1.82 \pm 0.01$ $+$ Kit Fast DNA MP $12.5 \pm 7.24^{a}$ $1.61 \pm 0.01$ $-$ PCI/SDS $1015.6 \pm 4.95^{b}$ $2.01 \pm 0.005$ $+$ PCI/SDS/NaCl $130.86 \pm 7,58^{\circ}$ $2.02 \pm 0.015$ $-$ Cooked shrimpCTAB $450.26 \pm 12.46^{\circ}$ $2.01 \pm 0.02$ $-$ CTAB $450.26 \pm 12.46^{\circ}$ $2.01 \pm 0.02$ $-$ TRITON $328.1 \pm 11.96^{\circ}$ $2.00 \pm 0.02$ $+$ DNeasy Tissue kit $87.17 \pm 3^{t}$ $1.71 \pm 0.01$ $+$ Kit Fast DNA MP $10.8 \pm 0.16^{\circ}$ $1.69 \pm 0.01$ $-$ PCI/SDS $1090.03 \pm 11.05^{b}$ $1.82 \pm 0$ $+$ PC/SDS $1090.03 \pm 11.05^{b}$ $1.82 \pm 0$ $+$ PCI/SDS/NaCl $59.51 \pm 7.39^{\circ}$ $1.73 \pm 0.01$ $-$ CTAB $62.3 \pm 1.20^{\circ}$ $1.97 \pm 0.03$ $+$ TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ $+$ DNeasy Tissue kit $64.13 \pm 8^{\circ}$ $1.41 \pm 0.01$ $+$ CTAB $62.3 \pm 1.20^{\circ}$ $1.71 \pm 0.005$ $-$ CTAB $62.3 \pm 1.20^{\circ}$ $1.71 \pm 0.005$ $-$ CTAB $62.3 \pm 1.20^{\circ}$ $1.92 \pm 0.015$ $+$ DNeasy Tissue kit $64.13 \pm 8^{\circ}$		PCI/SDS	$1679.5 \pm 26.52^{b}$	$1.98 \pm 0.015$	+
Plesionika heterocarpus         CHELEX $89.0 \pm 2.67^d$ $1.39 \pm 0.02$ -           CTAB $412.3 \pm 13.57^\circ$ $1.99 \pm 0.005$ +           TRITON $588.61 \pm 29.98^t$ $1.84 \pm 0.025$ +           DNeasy Tissue kit $84.17 \pm 0.6^d$ $1.82 \pm 0.01$ +           Kit Fast DNA MP $12.5 \pm 7.24^\circ$ $1.61 \pm 0.01$ -           PCI/SDS $1015.6 \pm 4.95^b$ $2.01 \pm 0.005$ +           PCI/SDS/NaCl $130.86 \pm 7,58^\circ$ $2.02 \pm 0.015$ -           CHELEX $53.03 \pm 11.16^4$ $1.13 \pm 0.02$ -           CTAB $450.26 \pm 12.46^\circ$ $2.01 \pm 0.02$ -           TRITON $328.1 \pm 11.96^\circ$ $2.00 \pm 0.02$ +           DNeasy Tissue kit $87.17 \pm 3^t$ $1.71 \pm 0.01$ +           Kit Fast DNA MP $10.8 \pm 0.16^\circ$ $1.69 \pm 0.01$ -           PCI/SDS $1090.03 \pm 11.05^b$ $1.82 \pm 0$ +           PCI/SDS/NaCl $59.51 \pm 7.39^\circ$ $1.73 \pm 0.01$ -           CTAB $62.3 \pm 1.20^\circ$ $1.97 \pm 0.03$ +           TRITON <td< td=""><td>PCI/SDS/NaCl</td><td><math>233.2 \pm 17.34^{\circ}</math></td><td><math>1.95 \pm 0.015</math></td><td>+</td></td<>		PCI/SDS/NaCl	$233.2 \pm 17.34^{\circ}$	$1.95 \pm 0.015$	+
$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$		CHELEX	$89.0 \pm 2.67^{d}$	$1.39 \pm 0.02$	_
		СТАВ	$412.3 \pm 13.57^{\circ}$	$1.99 \pm 0.005$	+
$ \begin{array}{ c c c c c c c c c c c c c c c c c c c$		TRITON	$588.61 \pm 29.98^{\rm f}$	$1.84 \pm 0.025$	+
$ \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$		DNeasy Tissue kit	$84.17 \pm 0.6$ <sup>d</sup>	$1.82 \pm 0.01$	+
$ \begin{tabular}{ c c c c c c c } \hline PCI/SDS & 1015.6 \pm 4.95^b & 2.01 \pm 0.005 & + \\ \hline PCI/SDS/NaCl & 130.86 \pm 7,58 ° & 2.02 \pm 0.015 & - \\ \hline CC1AB & 450.26 \pm 12.46° & 2.01 \pm 0.02 & - \\ \hline CTAB & 450.26 \pm 12.46° & 2.01 \pm 0.02 & - \\ \hline CTAB & 328.1 \pm 11.96° & 2.00 \pm 0.02 & + \\ \hline DNeasy Tissue kit & 87.17 \pm 3 f & 1.71 \pm 0.01 & + \\ \hline DNeasy Tissue kit & 87.17 \pm 3 f & 1.71 \pm 0.01 & + \\ \hline PCI/SDS & 1090.03 \pm 11.05^b & 1.69 \pm 0.01 & - \\ \hline PCI/SDS/NaCl & 59.51 \pm 7.39° & 1.73 \pm 0.01 & - \\ \hline PCI/SDS/NaCl & 59.51 \pm 7.39° & 1.73 \pm 0.01 & - \\ \hline CTAB & 62.3 \pm 1.20° & 1.97 \pm 0.03 & + \\ \hline TRITON & 391.9 \pm 17.8^d & 1.92 \pm 0.015 & + \\ \hline DNeasy Tissue kit & 64.13 \pm 8 ° & 1.41 \pm 0.01 & + \\ \hline Kit Fast DNA MP & 19.4 \pm 0.3^a & 1.71 \pm 0.005 & - \\ \hline PCI/SDS/NaCl & 47.35 \pm 8.70 ° & 1.36 \pm 0.005 & - \\ \hline PCI/SDS/NaCl & 47.35 \pm 8.70 ° & 1.36 \pm 0.005 & - \\ \hline PCI/SDS/NaCl & 47.35 \pm 8.70 ° & 1.36 \pm 0.005 & - \\ \hline CTAB & 52.7 \pm 37.83° & 1.50 \pm 0.015 & - \\ \hline CTAB & 52.7 \pm 37.83° & 1.50 \pm 0.01 $	Cooked shrimp	Kit Fast DNA MP	$12.5 \pm 7.24^{a}$	$1.61 \pm 0.01$	-
		PCI/SDS	$1015.6 \pm 4.95^{b}$	$2.01 \pm 0.005$	+
Cooked shrimpCHELEX $53.03 \pm 11.16^4$ $1.13 \pm 0.02$ $-$ CTAB $450.26 \pm 12.46^c$ $2.01 \pm 0.02$ $-$ TRITON $328.1 \pm 11.96^c$ $2.00 \pm 0.02$ $+$ DNeasy Tissue kit $87.17 \pm 3^{-1}$ $1.71 \pm 0.01$ $+$ Kit Fast DNA MP $10.8 \pm 0.16^a$ $1.69 \pm 0.01$ $-$ PCI/SDS $1090.03 \pm 11.05^b$ $1.82 \pm 0$ $+$ PCI/SDS/NaCl $59.51 \pm 7.39^c$ $1.73 \pm 0.01$ $-$ CTAB $62.3 \pm 12.0^c$ $1.97 \pm 0.03$ $+$ CTAB $62.3 \pm 12.0^c$ $1.97 \pm 0.03$ $+$ DNeasy Tissue kit $64.13 \pm 8^c$ $1.41 \pm 0.01$ $+$ DNeasy Tissue kit $64.13 \pm 8^c$ $1.41 \pm 0.01$ $+$ PCI/SDS $1079.53 \pm 7.16^{-b}$ $2.01 \pm 0.01$ $+$ PCI/SDS/NaCl $47.35 \pm 8.70^c$ $1.36 \pm 0.005$ $-$ PCI/SDS/NaCl $47.35 \pm 8.70^c$ $1.36 \pm 0.005$ $-$ CTAB $52.7 \pm 37.83^c$ $1.50 \pm 0.01$ $+$ PCI/SDS/NaCl $47.35 \pm 8.70^c$ $1.63 \pm 0.015$ $-$ DNeasy Tissue kit $58.87 \pm 6.2^c$ $1.52 \pm 0.01$ $+$		PCI/SDS/NaCl	$130.86 \pm 7,58$ °	$2.02 \pm 0.015$	-
$ \begin{array}{c c c c c c c c c c c c c c c c c c c $		CHELEX	$53.03 \pm 11.16^{d}$	$1.13 \pm 0.02$	-
$ \begin{array}{ c c c c c c c c c c c c c c c c c c c$		CTAB	$450.26 \pm 12.46^{\circ}$	$2.01 \pm 0.02$	-
DNeasy Tissue kit $87.17 \pm 3^{f}$ $1.71 \pm 0.01$ +Kit Fast DNA MP $10.8 \pm 0.16^{a}$ $1.69 \pm 0.01$ -PCI/SDS $1090.03 \pm 11.05^{b}$ $1.82 \pm 0$ +PCI/SDS/NaCl $59.51 \pm 7.39^{c}$ $1.73 \pm 0.01$ -CHELEX $44.53 \pm 15.10^{c}$ $1.21 \pm 0.01$ -CTAB $62.3 \pm 1.20^{c}$ $1.97 \pm 0.03$ +TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ +DNeasy Tissue kit $64.13 \pm 8^{c}$ $1.41 \pm 0.01$ +Kit Fast DNA MP $19.4 \pm 0.3^{a}$ $1.71 \pm 0.005$ -PCI/SDS $1079.53 \pm 7.16^{b}$ $2.01 \pm 0.01$ +PCI/SDS/NaCl $47.35 \pm 8.70^{c}$ $1.36 \pm 0.005$ -PCI/SDS/NaCl $47.35 \pm 8.70^{c}$ $1.36 \pm 0.005$ -CTAB $52,7 \pm 37.83^{c}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{c}$ $1.52 \pm 0.01$ +		TRITON	$328.1 \pm 11.96^{\circ}$	$2.00 \pm 0.02$	+
$ \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$		DNeasy Tissue kit	$87.17 \pm 3^{\text{f}}$	$1.71 \pm 0.01$	+
$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$	Frozen peeled precooked shrimp	Kit Fast DNA MP	$10.8 \pm 0.16^{a}$	$1.69 \pm 0.01$	-
Prozen peeled precooked shrimpPCI/SDS/NaCl $59.51 \pm 7.39^{\circ}$ $1.73 \pm 0.01$ $-$ CHELEX $44.53 \pm 15.10^{\circ}$ $1.21 \pm 0.01$ $-$ CTAB $62.3 \pm 1.20^{\circ}$ $1.97 \pm 0.03$ $+$ TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ $+$ DNeasy Tissue kit $64.13 \pm 8^{\circ}$ $1.41 \pm 0.01$ $+$ Kit Fast DNA MP $19.4 \pm 0.3^{a}$ $1.71 \pm 0.005$ $-$ PCI/SDS $1079.53 \pm 7.16^{b}$ $2.01 \pm 0.01$ $+$ PCI/SDS/NaCl $47.35 \pm 8.70^{\circ}$ $1.36 \pm 0.005$ $-$ CTAB $52.7 \pm 37.83^{\circ}$ $1.50 \pm 0.01$ $-$ TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ $-$ DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ $+$		PCI/SDS	$1090.03 \pm 11.05^{b}$	$1.82 \pm 0$	+
Frozen peeled precooked shrimpCHELEX $44.53 \pm 15.10^{\circ}$ $1.21 \pm 0.01$ -CTAB $62.3 \pm 1.20^{\circ}$ $1.97 \pm 0.03$ +TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ +DNeasy Tissue kit $64.13 \pm 8^{\circ}$ $1.41 \pm 0.01$ +Kit Fast DNA MP $19.4 \pm 0.3^{a}$ $1.71 \pm 0.005$ -PCI/SDS $1079.53 \pm 7.16^{b}$ $2.01 \pm 0.01$ +PCI/SDS/NaCl $47.35 \pm 8.70^{\circ}$ $1.36 \pm 0.005$ -CHELEX $91.3 \pm 19.52^{d}$ $1.14 \pm 0.015$ -CTAB $52.7 \pm 37.83^{\circ}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ +		PCI/SDS/NaCl	59.51 ± 7.39 <sup>e</sup>	$1.73 \pm 0.01$	-
CTAB $62.3 \pm 1.20^{\circ}$ $1.97 \pm 0.03$ +TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ +DNeasy Tissue kit $64.13 \pm 8^{\circ}$ $1.41 \pm 0.01$ +Kit Fast DNA MP $19.4 \pm 0.3^{\circ}$ $1.71 \pm 0.005$ -PCI/SDS $1079.53 \pm 7.16^{\circ}$ $2.01 \pm 0.01$ +PCI/SDS/NaCl $47.35 \pm 8.70^{\circ}$ $1.36 \pm 0.005$ -CHELEX $91.3 \pm 19.52^{d}$ $1.14 \pm 0.015$ -CTAB $52,7 \pm 37.83^{\circ}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ +		CHELEX	44.53 ± 15.10 °	$1.21 \pm 0.01$	-
TRITON $391.9 \pm 17.8^{d}$ $1.92 \pm 0.015$ +DNeasy Tissue kit $64.13 \pm 8^{\circ}$ $1.41 \pm 0.01$ +Kit Fast DNA MP $19.4 \pm 0.3^{a}$ $1.71 \pm 0.005$ -PCI/SDS $1079.53 \pm 7.16^{b}$ $2.01 \pm 0.01$ +PCI/SDS/NaCl $47.35 \pm 8.70^{\circ}$ $1.36 \pm 0.005$ -CHELEX $91.3 \pm 19.52^{d}$ $1.14 \pm 0.015$ -CTAB $52,7 \pm 37.83^{\circ}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ +		CTAB	62.3 ± 1.20°	$1.97 \pm 0.03$	+
DNeasy Tissue kit $64.13 \pm 8^{\circ}$ $1.41 \pm 0.01$ +Kit Fast DNA MP $19.4 \pm 0.3^{a}$ $1.71 \pm 0.005$ -PCI/SDS $1079.53 \pm 7.16^{b}$ $2.01 \pm 0.01$ +PCI/SDS/NaCl $47.35 \pm 8.70^{\circ}$ $1.36 \pm 0.005$ -CHELEX $91.3 \pm 19.52^{d}$ $1.14 \pm 0.015$ -CTAB $52,7 \pm 37.83^{\circ}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ +		TRITON	$391.9 \pm 17.8^{d}$	$1.92 \pm 0.015$	+
Kit Fast DNA MP $19.4 \pm 0.3^{a}$ $1.71 \pm 0.005$ -PCI/SDS $1079.53 \pm 7.16^{b}$ $2.01 \pm 0.01$ +PCI/SDS/NaCl $47.35 \pm 8.70^{c}$ $1.36 \pm 0.005$ -CHELEX $91.3 \pm 19.52^{d}$ $1.14 \pm 0.015$ -CTAB $52,7 \pm 37.83^{c}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{e}$ $1.52 \pm 0.01$ +		DNeasy Tissue kit	64.13 ± 8 °	$1.41 \pm 0.01$	+
PCI/SDS $1079.53 \pm 7.16^{b}$ $2.01 \pm 0.01$ +PCI/SDS/NaCl $47.35 \pm 8.70^{c}$ $1.36 \pm 0.005$ -CHELEX $91.3 \pm 19.52^{d}$ $1.14 \pm 0.015$ -CTAB $52,7 \pm 37.83^{c}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{c}$ $1.52 \pm 0.01$ +	Breaded Shrimp	Kit Fast DNA MP	$19.4 \pm 0.3^{a}$	$1.71 \pm 0.005$	-
PCI/SDS/NaCl $47.35 \pm 8.70^{\circ}$ $1.36 \pm 0.005$ -Breaded ShrimpCHELEX $91.3 \pm 19.52^{d}$ $1.14 \pm 0.015$ -CTAB $52,7 \pm 37.83^{\circ}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ +		PCI/SDS	1079.53 ± 7.16 <sup>b</sup>	$2.01 \pm 0.01$	+
Breaded ShrimpCHELEX $91.3 \pm 19.52^{d}$ $1.14 \pm 0.015$ -CTAB $52,7 \pm 37.83^{e}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{e}$ $1.52 \pm 0.01$ +		PCI/SDS/NaCl	47.35 ± 8.70 °	$1.36 \pm 0.005$	-
CTAB $52,7 \pm 37.83^{\circ}$ $1.50 \pm 0.01$ -TRITON $263.1 \pm 12.52^{f}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ +		CHELEX	$91.3 \pm 19.52^{d}$	$1.14 \pm 0.015$	-
TRITON $263.1 \pm 12.52^{\text{f}}$ $1.63 \pm 0.015$ -DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ +		CTAB	52,7 ± 37.83°	$1.50 \pm 0.01$	-
DNeasy Tissue kit $58.87 \pm 6.2^{\circ}$ $1.52 \pm 0.01$ +		TRITON	$263.1 \pm 12.52^{\text{f}}$	$1.63 \pm 0.015$	-
		DNeasy Tissue kit	58.87 ± 6.2 °	$1.52 \pm 0.01$	+



Journal abbreviation: J FisheriesSciences.com

**Figure 1:** PCR products from shrimp obtained with the UnivPF/UnivPR primers (315 bp) using different extraction methods. From left to right, Mq, molecular weight marker (100 pb DNA Ladder), DNA fragment products obtained with *Parapenaeus longirostris* (lane 1), cooked shrimp (lane 2), Frozen peeled precooked shrimp (lane3), *Metapenaeus monoceros* (lane 4), *Plezionika heterocarpus* (lane 5), *Metapenaeus kerathurus* (lane 6), Breaded shrimp (lane 7).



**Figure 2:** Loading plots from partial least square PLS1 regression carried out on all X-variables DO: A260/A280, DNA YIELD, EX1:Kit fast DNA method, EX2: PCI/SDS method, EX3: SDS/PCI/NaCl method, EX4: CHELEX method, EX5: CTAB/CI method, EX6: TRITON/CI method, EX7: DNeasy Tissue kit method and all shrimp sample foodstuff and the Y-variable (PCR RESULT: PCR amplification). All data in the PCA projections represent the means value of 6 replicates for each parameter.

#### Journal abbreviation: J FisheriesSciences.com

Similarly, both CTAB and Triton methods revealed considerable low DNA quantities, especially in breaded shrimp with poor DNA quality (A260/A280=1.50-1.63 respectively). In addition, poor reproducibility of DNA extraction was noticed as generally high standard errors were found for DNA quantification for fresh and processed shrimp. Among all trials, Fast DNA kit method yielded the lowest DNA recovery in all samples with values ranging from 10.8 to 41.7 ng DNA/µL of extract. Such kit was found inappropriate to DNA precipitation and recovery in shrimp matrix. The quantity and quality of the extracted DNA from shrimp clearly depended on the type of matrix and on the level of processing. Several parameters including the degree of DNA damage occurring during processing and the extraction procedure, as well as the average length of DNA fragments, affect considerably the success of DNA amplification. In this study, the DNA extracted from the different shrimp matrices were differently amplified with a set of primers amplifying the 315 bpcytochrome b fragments as shown in Table 1 and Figure 1. The amplification was successful particularly in fresh shrimp which was used as reference, but was difficult in breaded and cooked shrimp; proving that processing may damage the DNA. Thus DNA-fragments from fresh shrimp tissues were amplified by PCR independently of the extracting methods (Table 1). Concerning the DNA from the various processed shrimp matrices, only DNA extracted with the PCI/SDS and DNeasy Tissue kit gave a PCR product of the 315 bp-expected size for all the analysed samples (Figure 1 and Table 1).

#### Partial least square regression analysis

In this study, Partial least squares regression models (PLSR) were performed to identify the most significant variable relative to the different extraction method to PCR amplification of fresh and shrimp products (Martens & Næs, 1989). The PCR amplification was used as Y-matrix; all the others parameters (DO, A260/A280, DNA Yield, extraction methods and shrimp sample type) as X-variables.

The loading plot of shrimp (Figure 2) corresponded to a variance in PCR result of 62%, of which 59% was along PC1 and 3% along PC2. Of the total variation of all other variables (28%), 16% was along PC1 and 12% along PC2. Such result highlighted that fresh tissue (right upper region) is a discriminating parameter. The comparison of the different shrimp products revealed that, DO A260/A280, DNA YIELD, fresh shrimp (right hand side) and Breaded/Cooked shrimp (left hand side) contribute the most to important variations on PC1. The PC2 was mainly describing processed shrimp product. When comparing the extraction methods, results showed that DO, A260/A280, DNA YIELD, EX2, EX7, EX5 and EX6 (right hand side) and EX1, EX3, EX4 (left hand), contribute the most to the significant variation on PC1. Results showed that EX2, EX7, EX5 and EX6 were positively correlated with PCR amplification, while EX1, EX3, EX4 were negatively correlated with results of PCR. Such data confirmed results previously recorded for the DNA yield/purity and amplicon observed on agarose gel. In conclusion, comparing all products, the fresh shrimp yielded the highest amounts of DNA independently of the extraction methods; while matrices with high degree of processing such as breaded shrimp and cooked shrimp, yielded the lowest amounts of DNA.

The quality of DNA, with the highest amplicon length obtained from the fresh and the different shrimp food product, was best preserved using the SDS/PCI precipitation method, with the advantage of being the fastest and simplest protocol and is less expensive than the commercial kit. The relatively high values of the PLS1 of all the data (62%), suggest substantial difference between the extraction procedures and matrices when all the studied parameters were compared. Such statistical data confirm results previously recorded. Thus, the SDS/PCI method can be a good choice for processed shrimp products. Based on the present result, amplification of cyt b gene of a species particularly in processed food; dependent on the applied processing technique. When all studied parameters were compared, Partial least squares regression models (PLSR) suggested substantial difference between the extraction procedures and matrices; confirming previously recorded results. Among all tested protocols, the PCI/ SDS was identified as the most suitable method for the extraction of DNA from shrimp: with high DNA yields, a relative ease of use, less expensive than the commercial kit and the applicability for high throughput extraction from multiple specimens especially from highly processed foodstuffs.

#### **Conflict of Interests**

The authors declare no conflict of interest.

#### Acknowledgments

This work was conducted within the project entitled Biotechnologie Marine Vecteur d'Innovation et de Qualité - BIOVecQ PS1.3\_08 co-financed by the European Union through the European Neighborhood and Partnership Instrument (ENPI) Italy-Tunisia program and the Tunisian Ministry of Higher Education.

#### References

- Bauer, T., Weller, P., Hammes, W.P., Hertel, C. (2003) The effect of processing parameters on DNA degradation in food. Eur Food Research Tech 217, 338-343.
- Besbes, N., Fattouch, S., Sadok, S. (2011) Comparison of methods in the recovery and amplificability of DNA from fresh and processed sardine and anchovy muscle tissues. Food Chem 129, 665-671.
- Besbes, N, Jérôme, M., Bergé, J.P., Amri, M., Sadok, S. (2015) Authenticating the origin of different shrimp products on the Tunisian markets by PCR/RFLP method. African Journal of Biotechnology 14.
- Brown, W.M., Jr G.M., Wilson, A.C. (1979) Genetics Rapid evolution of animal mitochondrial DNA. Proc. Natl Acad Sci **76**, pp: 1967-1971.

- Cawthorn, D.M., Steinman, H.A., Witthuhn, R.C. (2011) Comparative study of different methods for the extraction of DNA from fish species commercially available in South Africa. Food Control **22**, pp: 231-244.
- Chapela, M.J., Sotelo, C.G., Ricardo, I.P.M., Miguel, A.P., Begoña, P.V., et al. (2007) Comparison of DNA extraction methods from muscle of canned tuna for species identification. Food Control 18, pp: 1211-1215.
- Chandrika, M., Maimunah, M., Zainon, M.N., Son, R. (2010) Identification of the species origin of commercially available processed food products by mitochondrial DNA analysis. International Food Research 17, pp: 867-876.
- Di Pinto, A., Forte, V., Guastadisegni, M.C., Martino, C., Schena, F.P. (2007) A comparison of DNA extraction methods for food analysis. Food Control 18, pp: 76-80.
- Hsieh, H.S., Chai, T.J., Hwang, D.F. (2005) Rapid PCR–RFLP method for the identification of 5 billfish species. J Food Science **70**, pp: 246-249.
- Jérôme, M., Lemaire, C., Bautista J.M., Fleurence, J.L., Etienne, M. (2003) Molecular Phylogeny and Species Identification of Sardines. J Agric Food Chem 51, pp: 7326-7332.
- Jérôme, M., Lemaire, C., Verrez-Bagnis, V., Etienne, M. (2003) Direct sequencing method for species identification of canned sardine and sardine-type products. J Agric Food Chem **51**, pp: 7326-7332.
- Kaniak-Golik A., Skoneczna A. (2015) Mitochondria–nucleus network for genome stability. Free Radical Biology and Medicine **2**, pp: 73-104.
- Lockley, A.K., Bardsley, R.G. (2000) DNA-based methods for food authentication. Trend Food Sci 11, pp: 67-77.
- Martens, H., Næs, T. (1989) Multivariate calibration. Chichester, UK: Wiley.

- Olexová, L., Dovicovicová, L., Kuchta, T. (2004). Comparison of three types of methods for the isolation of DNA from flours, biscuits and instant paps. Eur Food Research Tech **218**, pp: 390-393.
- Pascoal, A., Barros-Velazquez, J., Cepeda, A., Gallardo, J.M., Calo-Mata, P. (2008) A polymerase chain reaction-restriction fragment length polymorphism method based on the analysis of a 16S rRNA/tRNA(Val) mitochondrial region for species identification of commercial penaeid shrimps (Crustacea: Decapoda: Penaeoidea) of food interest. Electrophoresis 29, pp: 499-509.
- Prado, M., Pascoal, A., Calo-Mata, P., Barros-Velázquez, J. (2005). Molecular methods for meat species identification and its application to foodstuffs and animal feeds. In A. P. Riley (Ed.), Food research, safety and policies. New York: Nova Science Publishers Inc pp: 199-221.
- Quinta, R., Gomes, L., Teia dos Santos, A. (2004) A mitochondrial DNA PCR-RFLP marker for population studies of the black scabbard fish (*Aphanopus carbo*). J. Marine Sci **61**, pp: 864-867.
- Rumpho, M.E., Mujer, C.V., Andrews, D.L., Manhart, J.R., Pierce, S.K. (1994) Extraction of DNA from mucilaginous tissues of a sea slug (*Elysia chlorotica*). Biotechniques 17, pp: 1097-1101.
- Sambrock, J., Fritsch, E.F., Maniatis, T. (1989) Molecular cloning: a laboratory manual (2nd ed.). Cold Spring Harbor, New York: Cold Spring Harbor Laboratory press.
- Wiesner, R.J., Swift, H., Zak, R. (1991) Purification of mitochondrial DNA from total cellular DNA of small tissue samples. Gene. **98**, pp: 277-281.